

Enhanced *Laboratory Portfolio*



Fisherbrand glassware quality ensured with:

- High clarity, uniform thickness and durable glass
- Long lasting graduation marking
- Dedicated QA for calibration and testing
- Superior packaging for safe handling

Fisherbrand®

Value + Added

- Easily access Certification of Conformity (COC) for FisherBrand® Class “A” glassware via online portal (<https://portal.thermofisher.co.in/CoC/viewCoc.aspx>) after filling necessary details like Product code & Unit no.




- *Flasks*
- *Cylinders*
- *Pipettes*
- *Burettes*

Note: COC for "A" class Volumetric ware that have reached our warehouse after 1st August, 2015, will be available on portal. For Volumetric ware received before 1st August, 2015, online portal allows customer to send request for soft copy on their email-id.




Volumetric Glassware

- Precision Bore Tubing
- Amber Graduation
- Individually Calibrated
- Complete Traceability
- ISO & ASTM specifications
- NABL Certified




Aspirator Bottle with GL45 and tube at Bottom




Part No.	Capacity (ml)	App. O.D x Height (mm)	Neck Size	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6457A	250	70 x 138	GL 45	6	341.33	2048
FBG6457B	500	86 x 176	GL 45	6	385.33	2312
FBG6457C	1000	101 X 225	GL 45	4	694.50	2778
FBG6457D	2000	136 X 260	GL 45	2	1090.50	2181

Aspirator Bottle with GL 45 and socket for Stopper



Part No.	Capacity (ml)	App. O.D x Height (mm)	Neck Size	Outlet Ground Joint Size	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6458A	1000	101 X 225	GL 45	24/29	4	907.75	3631
FBG6458B	2000	136 X 260	GL 45	29/32	2	1277.00	2554
FBG6458C	5000	180 X 330	GL 45	29/32	1	7522.00	7522
FBG6458D	10000	227 X 410	GL 45	29/32	1	11139.00	11139

Aspirator Bottle with GL 45 and Stopcock




Part No.	Capacity (ml)	App. O.D x Height (mm)	Neck Size	Outlet Ground Joint Size	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6459A	1000	101 X 225	GL 45	24/29	4	2534.75	10139
FBG6459B	2000	136 X 260	GL 45	29/32	2	2935.00	5870
FBG6459C	5000	181 X 330	GL 45	29/32	1	9895.00	9895
FBG6459D	10000	227 X 410	GL 45	29/32	1	13419.00	13419

Burette With Straight Bore Moulded Screw Thread Glass Stopcock, Class A NABL Certified




Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63055N	10	0.05	0.02	1	4035.00	4035
FBG63056N	25	0.10	0.05	1	4191.00	4191
FBG63057N	50	0.10	0.05	1	4786.00	4786
FBG63058N	100	0.20	0.10	1	8674.00	8674

Burette With Straight Bore Ptfе Key Stopcock, Class A NABL Certified



Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63075N	10	0.05	0.02	1	3917.00	3917
FBG63076N	25	0.10	0.05	1	4047.00	4047
FBG63077N	50	0.10	0.05	1	4584.00	4584
FBG63078N	100	0.20	0.10	1	6705.00	6705

Measuring Cylinder, With Hexagonal Glass Base, Class A NABL Certified




Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62593N	5	0.1	0.05	115	12	1082.17	12986
FBG62594N	10	0.2	0.10	140	2	1196.50	2393
FBG62596N	25	0.5	0.25	170	2	1402.50	2805
FBG62597N	50	1.0	0.50	200	2	1668.50	3337
FBG62598N	100	1.0	0.50	260	2	1894.00	3788
FBG6259BN	250	2.0	1.00	300	2	2872.00	5744
FBG6259CN	500	2.0	2.50	390	2	3919.00	7838
FBG6259EN	1000	5.0	5.00	470	1	4705.00	4705
FBG6259FN	2000	10.0	10.00	620	1	9081.00	9081



Glassware


Flask Volumetric, With Interchangeable Glass Stopper, Class A NABL Certified



Part No.	Capacity (ml)	Tolerance (\pm ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62194N	5	0.025	10/19	6	1973.83	11843
FBG62195N	10	0.025	10/19	6	1706.83	10241
FBG62192N	20	0.04	10/19	2	1634.50	3269
FBG62196N	25	0.04	10/19	2	1590.50	3181
FBG62197N	50	0.06	12/21	2	1590.50	3181
FBG62198N	100	0.10	14/23	2	1548.00	3096
FBG6219AN	200	0.15	14/23	2	2005.00	4010
FBG6219BN	250	0.15	14/23	2	2025.00	4050
FBG6219CN	500	0.25	19/26	2	2689.50	5379
FBG6219EN	1000	0.40	24/29	2	3904.50	7809
FBG6219FN	2000	0.60	29/32	1	6323.00	6323


* Flask in Test Tube form

Pipette Graduated, Mohr Type, Color Coded, Class A NABL Certified




Part No.	Capacity (ml)	Sub division (ml)	Tolerance (\pm ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63341N	1	0.01	0.006	1	836.00	836
FBG63342N	2	0.02	0.010	1	876.00	876
FBG63343N	5	0.05	0.030	1	1035.00	1035
FBG63344N	10	0.10	0.050	1	1095.00	1095
FBG63346N	25	0.20	0.100	1	1278.00	1278

Pipette Volumetric With One Mark, Bulb Type, Class A NABL Certified




Part No.	Capacity (ml)	Tolerance (\pm ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63321N	1	0.007	1	607.00	607
FBG63322N	2	0.010	1	627.00	627
FBG63323N	5	0.015	1	679.00	679
FBG63324N	10	0.020	1	834.00	834
FBG63325N	20	0.030	1	1089.00	1089
FBG63326N	25	0.030	1	1174.00	1174
FBG63327N	50	0.050	1	1310.00	1310
FBG63328N	100	0.080	1	2055.00	2055

Weighing Scoop



Part No.	Capacity (ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG65048	1 ml	10	223.10	2231
FBG65049	2 ml	10	223.10	2231
FBG65050	3 ml	10	195.40	1954
FBG65051	6 ml	10	223.10	2231
FBG65052	10 ml	10	257.00	2570


Beakers, Griffin (Low Form), Borosilicate Glass, Graduated with Spout, in compliance with ISO



Part No.	Capacity (ml)	App. O.D (mm)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63003	5*	22	30	10	78.40	784
FBG63004	10*	26	35	10	67.00	670
FBG63006	25	34	50	10	78.40	784
FBG63007	50	42	60	10	66.00	660
FBG63008	100	50	70	10	70.00	700
FBG63009	150	60	80	10	70.00	700
FBG6100B	250	70	95	10	72.00	720
FBG6100C	500	83	115	8	116.00	928
FBG6100D	600	90	125	8	150.00	1200
FBG6100E	1000	105	145	6	243.00	1458
FBG6100F	2000	132	185	2	526.50	1053
FBG6100G	3000	146	214	1	1231.00	1231
FBG6100J	5000	170	265	1	1787.00	1787
FBG6100K	10000	220	320	1	6173.00	6173

* non graduated

Beaker, Berzelius (Tall Form), Graduated with Spout, in compliance with ISO



Part No.	Capacity (ml)	App. O.D (mm)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63027	50	38	70	10	108.30	1083
FBG63028	100	48	80	10	123.10	1231

Bottles, Reagent, Plain, Narrow Mouth with I/C Stopper, in compliance with ISO




Part No.	Capacity (ml)	Dimension App. O.D x Height	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62037	60	45 x 93	14/23	2	311.50	623
FBG62036	100	50 x 98	19/26	2	333.50	667
FBG62038	125	55 x 93	19/26	2	296.00	592
FBG6203B	250	64 x 147	19/26	2	330.00	660
FBG6203C	500	80 x 172	24/29	6	446.33	2678
FBG6203E	1000	104 x 200	29/32	4	610.75	2443
FBG6203F	2000	129 x 265	34/35	2	1033.00	2066

Bottles, Reagent, Amber, Narrow Mouth with I/C Stopper, in compliance with ISO



Part No.	Capacity (ml)	Dimension App. O.D X Height	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62057	60	45 x 93	14/23	2	451.00	902
FBG62056	100	50 x 98	19/26	2	471.00	942
FBG62058	125	55 x 93	19/26	2	471.00	942
FBG6205B	250	64 x 147	19/26	2	567.00	1134
FBG6205C	500	80 x 172	24/29	6	665.00	3990
FBG6205E	1000	104 x 200	29/32	4	1051.50	4206
FBG6205F	2000	129 x 265	34/35	2	1784.00	3568

Media Bottles, Borosilicate Glass with Screw Cap & Pouring Ring, Autoclavable



Part No.	Capacity (ml)	Dimension (Mm) App. O.D X Height	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FS6439525	25*	36 x 70	20	158.40	3168
FS6439550	50*	46 x 87	20	182.60	3652
FB800100	100	56 x 100	10	363.50	3635
FB800250	250	70 x 138	10	362.90	3629
FB800500	500	86 x 176	10	407.70	4077
FB8001000	1000	101 x 196	10	653.00	6530
FB8002000	2000	136 x 252	10	1123.80	11238
FB8003000	3000	158 x 295	1	5106.00	5106
FB8005000	5000	183 x 320	6	6915.67	41494
FB80010000	10000	235 x 400	1	15006.00	15006
FB33182	20000	296 x 495	1	10585.00	10585

* This bottle is without pouring ring.

Media Bottles, Borosilicate Glass, Amber Media with Screw Cap & Pouring Ring Autoclavable



Part No.	Capacity (ml)	Dimension (mm) App. O.D x Height	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FSG3546	25*	36 x 70	20	356.40	7128
FSG3547	50*	46 x 87	20	409.20	8184
FB800100A	100	56 x 100	10	567.80	5678
FB800250A	250	70 x 138	10	645.20	6452
FB800500A	500	86 x 176	10	880.80	8808
FB8001000A	1000	101 x 196	10	1175.10	11751
FB8002000A	2000	136 x 252	10	2009.50	20095
FB8003000A	3000	158 x 295	1	9880.00	9880
FB8005000A	5000	183 x 320	1	8086.00	8086
FB80010000A	10000	235 x 400	1	9346.00	9346
FB33182A	20000	296 x 495	1	21692.00	21692

* This bottle is without pouring ring.

Bottles, B.O.D, with I/C Stopper, specified by American Public Health Association



Part No.	Capacity (ml)	Dimension (mm) App. O.D x Height	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62558	125	55 x 120	19/26	2	434.00	868
FBG6255B	300	70 X 150	19/26	2	533.00	1066

Bottles, Relative Density, With Capillary Bore PTFE Stopper




Part No.	Capacity (ml)	Tolerance (± ml)	Dimensions App. O.D x Ht	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62565	10	1	24 x 46	10/15	10	415.50	4155
FBG62566	25	2	38 x 56	10/15	10	451.00	4510
FBG62567	50	3	49 x 66	10/15	10	482.90	4829



Glassware


Bottles, Weighing, with Ground in Stopper



Part No.	Capacity (ml)	App. O.D (mm)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63043	5	20	40	10	252.00	2520
FBG63044	15	25	60	10	266.00	2660
FBG63045	20*	50	35	10	553.00	5530
FBG63046	25	30	65	10	325.00	3250
FBG63042	40*	60	30	10	781.00	7810
FBG63047	60	40	90	10	494.00	4940


* Squat form

Bottles Solution, Tooled Neck, Plain



Part No.	Capacity (ml)	Dimension(mm) App. O.D x Height x Neck ID	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6257B	250	64 x 147 x 18	6	266.17	1597
FBG6257C	500	80 x 172 x 23	6	312.17	1873
FBG6257E	1000	104 x 200 x 28	4	511.50	2046
FBG6257F	2000	129 x 265 x 33	2	909.00	1818
FBG6257G	3000	158 x 260 x 44	1	2106.00	2106
FBG6257J	5000	190 x 320 x 44	1	5277.00	5277
FBG6257K	10000	215 x 410 x 54	1	9148.00	9148
FBG6257L	20000	300 x 485 x 54	1	14587.00	14587

Bottles, Solution, Tooled Neck, Amber




Part No.	Capacity (ml)	Dimension(mm) App. O.D x Height x Neck ID	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6357B	250	64 x 147 x 18	6	676.33	4058
FBG6357C	500	80 x 172 x 23	6	863.00	5178
FBG6357E	1000	104 x 200 x 28	4	1451.00	5804
FBG6357F	2000	129 x 265 x 33	2	1624.50	3249
FBG6357K	10000	215 x 410 x 54	1	16204.00	16204
FBG6357L	20000	300 x 485 x 54	1	20796.00	20796

Burettes, Straight Bore Stopcock, Glass Key, Class 'A', in compliance with ISO



Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63055	10	0.05	0.02	4	1770.75	7083
FBG63056	25	0.10	0.05	4	1770.00	7080
FBG63057	50	0.10	0.05	4	3891.00	15564
FBG63058	100	0.20	0.10	4	2948.75	11795

Burettes, Straight Bore Stopcock, Glass Key, Class 'B', Amber Graduation, in compliance with ISO



Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62065	10	0.05	0.05	4	631.00	2524
FBG62066	25	0.10	0.10	4	616.00	2464
FBG62067	50	0.10	0.10	4	655.00	2620
FBG62068	100	0.20	0.20	4	940.75	3763

Burettes, Rotaflow Screw Thread Stopcock With PTFE Key, Class 'A' Amber Stain Graduation, in compliance with ISO



Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63075	10	0.05	0.02	4	1534.50	6138
FBG63076	25	0.10	0.05	4	3309.00	13236
FBG63077	50	0.10	0.05	4	1922.00	7688
FBG63078	100	0.20	0.10	4	2675.00	10700

Burettes, Rotaflow Screw Thread Stopcock With PTFE Key, Class 'B' Amber Stain Graduation, in compliance with ISO



Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62075	10	0.05	0.05	4	560.75	2243
FBG62076	25	0.10	0.10	4	844.00	3376
FBG62077	50	0.10	0.10	4	932.00	3728
FBG62078	100	0.20	0.20	4	851.50	3406

Burettes, Straight Bore Stopcock PTFE Key, in compliance with ISO




Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62069	50	0.10	0.10	4	787.00	3148

Cylinders, Nessler's Flat Bottom, Graduated




Part No.	Capacity (ml)	App. O.D. (mm)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62087	50	26	150	10	160.00	1600
FBG62088	100	34	180	10	238.00	2380

Cylinders, Graduated, Hexagonal Base and Spout Class 'A', in compliance with ISO




Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62593	5	0.1	0.05	115	2	403.00	806
FBG62594	10	0.2	0.10	140	2	461.00	922
FBG62596	25	0.5	0.25	170	2	529.00	1058
FBG62597	50	1.0	0.50	200	2	640.00	1280
FBG62598	100	1.0	0.50	260	2	727.00	1454
FBG62599	150	2.0	1.00	300	2	966.00	1932
FBG6259B	250	2.0	1.00	335	2	1072.00	2144
FBG6259C	500	5.0	2.50	390	2	1494.00	2988
FBG6259E	1000	10.0	5.00	470	2	1829.00	3658
FBG6259F	2000	20.0	10.00	620	1	3556.00	3556

Cylinders, Graduated, Hexagonal Base And Spout, Class 'B', in compliance with ISO



Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62584	10	0.2	0.2	2	257.00	514
FBG62586	25	0.5	0.5	2	305.00	610
FBG62587	50	1.0	1.0	2	330.00	660
FBG62588	100	1.0	1.0	2	378.00	756
FBG6258B	250	2.0	2.0	1	699.00	699
FBG6258C	500	5.0	5.0	1	912.00	912
FBG6258E	1000	10.0	10.0	1	1383.00	1383

Cylinders, Graduated, Hexagonal Base and Spout, Class 'A', in compliance with ASTM (New)



Part No.	Capacity (ml)	Graduations (ml)	Tolerance (± ml)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63593	5	0.1	0.05	115	2	466.00	932
FBG63594	10	0.2	0.10	140	2	525.50	1051
FBG63596	25	0.5	0.17	170	2	605.50	1211
FBG63597	50	1.0	0.25	200	2	734.50	1469
FBG63598	100	1.0	0.50	260	2	808.00	1616
FBG6359B	250	2.0	1.00	335	2	1208.00	2416
FBG6359C	500	5.0	2.00	390	2	1591.00	3182
FBG6359E	1000	10.0	3.00	470	2	1960.00	3920
FBG6359F	2000	20.0	6.00	620	1	3677.00	3677




Glassware

Cylinders, Graduated, Round Base And Spout, Class 'B', in compliance with ISO




Part No.	Capacity (ml)	Graduations (ml)	Tolerance (\pm ml)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62094	10	0.2	0.2	140	2	247.00	494
FBG62096	25	0.5	0.5	170	2	303.00	606
FBG62097	50	1.0	1.0	200	2	325.00	650
FBG62098	100	1.0	1.0	260	2	414.50	829
FBG6209B	250	2.0	2.0	335	2	730.00	1460
FBG6209C	500	5.0	5.0	390	1	918.00	918
FBG6209E	1000	10.0	10.0	470	1	1408.00	1408
FBG6209F	2000	10.0	10.0	570	1	2406.00	2406

Cylinder, Crow, With I/C Stopper, in compliance with ISO



Part No.	Capacity (ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62807	50	2	447.50	895
FBG62808	100	2	536.50	1073

Cylinders, Penny Head with I/C Stopper, Amber Stain, Class 'A', in compliance with ISO




Part No.	Capacity (ml)	Graduations (ml)	Tolerance (\pm ml)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62504	10	0.2	0.1	160	2	733.00	1466
FBG62506	25	0.5	0.25	186	2	844.00	1688
FBG62507	50	1.0	0.5	206	2	955.00	1910
FBG62508	100	1.0	0.5	285	2	1121.00	2242
FBG6250B	250	2.0	1.0	370	1	1761.00	1761
FBG6250A	500	5.0	2.5	430	1	2261.00	2261
FBG6250E	1000	10.0	5.0	500	1	2770.00	2770

Extraction Apparatus, Soxhlet, with Allihn Condenser and I/C Joints




Part No.	Capacity (ml)	Flask size (mm)	Extractor socket (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62727	60	150	34/35	1	1348.00	1348
FBG62728	100	250	40/38	1	1573.00	1573
FBG6272A	200	500	50/42	1	2595.00	2595

Cylinders, Graduated, With I/C Stopper, Penny Head, Round Base and Spout, Class 'B', in compliance with ISO




Part No.	Capacity (ml)	Graduations (ml)	Tolerance (\pm ml)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62604	10	0.2	0.2	160	2	378.00	756
FBG62606	25	0.5	0.5	186	2	422.00	844
FBG62607	50	1.0	1.0	206	2	476.00	952
FBG62608	100	1.0	1.0	285	2	582.00	1164
FBG6260B	250	2.0	2.0	370	1	985.00	985
FBG6260C	500	5.0	5.0	430	1	1368.00	1368
FBG6260E	1000	10.0	10.0	500	1	1975.00	1975

Extractors, Soxhlet, with I/C Joints



Part No.	Capacity (ml)	I/C joint top (mm)	I/C joint bottom (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62737	60	34/35	24/29	2	724.50	1449
FBG62738	100	40/38	29/32	2	951.00	1902
FBG6273A	200	50/42	24/29	2	1114.00	2228

Condenser, Allihn, for Extraction Apparatus



Part No.	Capacity (ml)	I/C joint bottom (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62747	60	34/35	2	669.50	1339
FBG62748	100	40/38	2	691.50	1383
FBG6274A	200	50/42	2	946.50	1893

Condenser, Graham, with Drip Tip, I/C Joints



Part No.	Length (mm)	I/C joint Inner	I/C joint Outer	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62773	300	24/29	24/29	1	864.00	864
FBG62774	400	24/29	24/29	1	947.00	947
FBG62775	500	24/29	24/29	1	1031.00	1031

Condenser, Liebig, with Drip Tip, I/C Joints




Part No.	Length (mm)	I/C joint Inner	I/C joint Outer	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62752	200	19/26	19/26	1	355.00	355
FBG62753	300	24/29	24/29	1	428.00	428
FBG62754	400	24/29	24/29	1	473.00	473

Desiccator with Cover, Knob Top, made from 3.3 Borosilicate Glass




Part No.	Approx. I.D. of Ground Flange (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6279A	100	1	5887.00	5887
FBG6279B	150	1	6283.00	6283
FBG6279C	200	1	7101.00	7101
FBG6279D	250	1	10323.00	10323
FBG6279E	300	1	20911.00	20911

Desiccator with Cover, Knob Top, made from 3.3 Borosilicate Glass




Part No.	Approx. I.D. of Ground Flange (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6379A	100	1	6947.00	6947
FBG6379B	150	1	7387.00	7387
FBG6379C	200	1	9202.00	9202
FBG6379D	250	1	14245.00	14245
FBG6379E	300	1	23367.00	23367

Condenser, Allihn, with Drip Tip, I/C Joints



Part No.	Length (mm)	I/C joint Inner	I/C joint Outer	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62762	200	19/26	19/26	1	477.00	477
FBG62763	300	24/29	24/29	1	609.00	609
FBG62764	400	24/29	24/29	1	721.00	721
FBG62766	600	29/32	29/32	1	1305.00	1305

Desiccator, Vacuum, with PTFE Stopcock, made from 3.3 Borosilicate Glass




Part No.	Approx. I.D. of Ground Flange (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6289A	100	1	13957.00	13957
FBG6289B	150	1	20280.00	20280
FBG6289C	200	1	21339.00	21339
FBG6289D	250	1	35521.00	35521
FBG6289E	300	1	36309.00	36309




Glassware

Distilling Apparatus, Ammonia with Graham Condenser



Part No.	Capacity (ml)	Interchangeable Flask	Joint (mm) Condenser	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62781	500	24/29	19/26	1	2694.00	2694


Flasks, Boiling, Flat Bottom, in compliance with ISO



Part No.	Capacity (ml)	Approx. O.D. x Height (mm)	Approx. neck O.D. (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62126	25	41 x 90	20	10	88.20	882
FBG62127	50	51 x 100	26	10	88.90	889
FBG62128	100	64 x 110	26	10	109.80	1098
FBG62129	150	75 x 132	28	10	110.80	1108
FBG6212B	250	85 x 140	34	10	121.80	1218
FBG6212A*	300	87 x 160	24	10	194.30	1943
FBG6212C	500	105 x 171	33	6	313.83	1883
FBG6212E	1000	131 x 200	42	6	335.33	2012
FBG6212F	2000	166 x 250	50	2	1836.00	3672
FBG6212G	3000	186 x 300	50	1	3222.00	3222
FBG6212J	5000	223 x 317	60	1	2696.00	2696
FBG6212K	10000	279 x 400	60	1	6227.00	6227

* Suitable for use in Polenske Apparatus.

Flasks, Boiling, Round Bottom, in compliance with ISO




Part No.	Capacity (ml)	Approx. O.D. x Height (mm)	Approx. neck O.D. (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62116	25	41 x 90	20	10	86.90	869
FBG62117	50	51 x 95	26	10	88.90	889
FBG62118	100	64 x 110	26	10	109.80	1098
FBG62119	150	75 x 137	28	10	110.60	1106
FBG6211B	250	85 x 144	34	10	110.40	1104
FBG6211C	500	105 x 175	33	10	188.20	1882
FBG6211E	1000	131 x 200	42	6	323.00	1938
FBG6211F	2000	166 x 260	50	6	679.33	4076
FBG6211G	3000	190 x 310	50	1	1736.00	1736
FBG6211J	5000	225 x 350	60	1	2933.00	2933
FBG6211K	10000	279 x 400	65	1	5601.00	5601

Flasks, Boiling, Flat Bottom, with I/C Joint, in compliance with ISO



Part No.	Capacity (ml)	I/C joint (mm)	App. Height (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62629	150	24/29	115	10	166.20	1662
FBG6262B	250	24/29	130	10	177.00	1770
FBG6262C	500	24/29	148	10	266.40	2664

Flasks, Boiling, Round Bottom, I/C Joint, in compliance with ISO



Part No.	Capacity (ml)	I/C joint (mm)	App. Height (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62617	50	19/26	90	10	149.90	1499
FBG62618	100	24/29	110	10	155.30	1553
FBG62619	150	24/29	115	10	183.40	1834
FBG6261B	250	24/29	130	10	199.60	1996
FBG6261C	500	24/29	150	6	277.50	1665
FBG6261E	1000	24/29	188	2	389.50	779


Flasks, Boiling, R B, Two Necks, Centre & Angled Side Necks



Part No.	Capacity (ml)	Interchangeable Centre Socket	Joint (mm) Side Socket	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62638	100	24/29	19/26	1	244.00	244
FBG6263B	250	24/29	19/26	1	332.00	332
FBG6263C	500	24/29	19/26	1	456.00	456

* For any other neck size contact us.

Flasks, Boiling, R B, Two Necks, Centre & Parallel Side Necks



Part No.	Capacity (ml)	Interchangeable Centre Socket	Joint (mm) Side Socket	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62648	100	24/29	19/26	1	244.00	244
FBG6264B	250	24/29	19/26	1	332.00	332
FBG6264C	500	24/29	19/26	1	455.00	455

* For any other neck size contact us.


Flasks, Boiling, R B, Three Necks, Centre & Two Angled Side Necks



Part No.	Capacity (ml)	Interchangable Centre Socket	Joint (mm) Side Socket	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62658	100	24/29	19/26	1	305.00	305
FBG6265B	250	24/29	19/26	1	422.00	422
FBG6265C	500	24/29	19/26	1	642.00	642
FBG6265E	1000	24/29	19/26	1	779.00	779

* For any other neck size contact us.


Flasks, Boiling, R B, Three Necks, Centre & Two Parallel Side Necks



Part No.	Capacity (ml)	Interchangable Centre Socket	Joint (mm) Side Socket	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62668	100	24/29	19/26	1	333.00	333
FBG6266B	250	24/29	19/26	1	433.00	433
FBG6266C	500	24/29	19/26	1	612.00	612
FBG6266E	1000	24/29	19/26	1	779.00	779
FBG6266F	2000	24/29	19/26	1	1687.00	1687
FBG6266J	5000	24/29	19/26	1	5436.00	5436


* For any other neck size contact us.

Haffkine Flask (Culture Flask)




Part No.	Capacity (ml)	App. O.D x Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6218G	3000	241 x 273	1	2832.00	2832
FBG6218H	4000	245 x 310	2	3277.00	6554

Flask, Distilling, with Side Arm



Part No.	Capacity (ml)	App.O.D. x Height (mm)	Neck O.D. (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62698	100	65 x 215	21	1	212.00	212
FBG62699	125	68 x 214	19	1	234.00	234
FBG6269B	250	85 x 226	25	1	267.00	267
FBG6269C	500	105 x 240	28	1	358.00	358

Flask, Erlenmeyer (Conical Flask), with I/C Joint




Part No.	Capacity (ml)	App. Height (mm)	I/C Joint (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG6216A	100	90	24/25	10	199.60	1996
FBG6216B	250	145	24/29	10	210.50	2105
FBG6216C	500	175	24/29	10	299.30	2993

Flask, Erlenmeyer, Conical, with I/C Stopper



Part No.	Capacity (ml)	App. Height (mm)	I/C Joint (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62676	25	65	14/15	10	179.00	1790
FBG62677	50	85	19/20	10	213.00	2130
FBG62678	100	90	24/25	10	233.40	2334
FBG6267B	250	145	29/25	10	277.30	2773
FBG6267C	500	175	29/25	10	400.10	4001
FBG6267E	1000	175	29/25	6	558.00	3348

Flask, Conical, with Screw Cap




Part No.	Capacity (ml)	App. Height (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FS61138	100	100	20	269.75	5395
FS61139	150	125	20	318.15	6363
FS6113B	250	145	20	385.60	7712
FS6113C	500	180	10	465.90	4659
FS6113E	1000	225	10	867.20	8672




Glassware

Flask, Conical, Amber, with Screw Cap




Part No.	Capacity (ml)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FS62138	100	100	20	512.75	10255
FS62139	150	125	20	577.20	11544
FS6213B	250	145	20	661.10	13222
FS6213C	500	180	10	763.70	7637
FS6213E	1000	225	10	1990.50	19905

Flasks, Erlenmeyer, Conical, Graduated, in compliance with ISO




Part No.	Capacity (ml)	App. O.D X height (mm)	Neck O.D. (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG61134	10	30 x 55	15	10	69.10	691
FBG61136	25	42 x 70	16	10	77.90	779
FBG61137	50	51 x 82	20	10	77.00	770
FBG61138	100	64 x 105	23	10	77.00	770
FBG61139	150	72 x 124	25	10	97.00	970
FBG6113B	250	85 x 140	30	10	101.00	1010
FBG6113C	500	104 x 176	35	10	373.00	3730
FBG6113E	1000	131 x 220	40	6	175.00	1050
FBG6113F	2000	165 x 286	48	2	558.00	1116
FBG6113G	3000	165 x 286	48	1	1397.00	1397
FBG6113J	5000	220 x 377	50	1	2356.00	2356

Flasks, Erlenmeyer, Conical, Graduated, Amber, in compliance with ISO




Part No.	Capacity (ml)	App. O.D X height (mm)	Neck O.D. (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62158	100	64 x 105	23	10	150.00	1500
FBG62159	150	72 x 124	25	10	175.00	1750
FBG6215B	250	85 x 140	30	10	209.00	2090
FBG6215C	500	104 x 176	35	10	350.00	3500
FBG6215E	1000	131 x 220	40	10	592.00	5920

Flask, Filtration, Heavy Wall, with Outlet for Vacuum



Part No.	Capacity (ml)	Approx. O.D x Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6268C	500	105 x 184	2	511.00	1022
FBG6268E	1000	136 x 250	2	762.00	1524
FBG6268F	2000	168 x 295	2	1531.50	3063
FBG6268G	5000	237 x 450	2	6254.50	12509

Flasks, Iodine Determination, With Stopper



Part No.	Capacity (ml)	Joint mm (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6214A	50	24/29	2	197.50	395
FBG6214D	100	24/29	2	223.00	446
FBG6214B	250	24/29	2	252.00	504
FBG6214C	500	24/29	2	299.00	598

Flasks, Kjeldahl, Round Bottom, Long Neck



Part No.	Capacity (ml)	App. O.D (mm)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6217B	430	81	300	2	290.00	580
FBG6217C	520	97	330	2	592.50	1185
FBG6117D	800	113	375	1	693.00	693

Flasks, Kjeldahl, Round Bottom, I/C Joint And Long Neck




Part No.	Capacity (ml)	App. O.D (mm)	App. Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6218C	500	101	325	2	715.00	1430
FBG6118D	800	115	350	1	782.00	782

Flask, Volumetric, Sugar Estimation



Part No.	Capacity (ml)	Joint mm (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62218	100	0.15	2	313.00	626

Flasks, Volumetric, Class 'A', with Stopper, Made of 3.3 Borosilicate Glass, - USP, in compliance with ASTM




Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63190S	1*	0.01	10/19	2	874.00	1748
FBG63191S	2*	0.01	10/19	2	874.00	1748
FBG63193S	3	0.01	10/19	2	874.00	1748
FBG63194S	5	0.02	10/19	2	723.00	1446
FBG63195S	10	0.02	10/19	2	704.00	1408
FBG63192S	20	0.03	10/19	2	703.00	1406
FBG63196S	25	0.03	10/19	2	650.00	1300
FBG63197S	50	0.05	10/19	2	621.00	1242
FBG63198S	100	0.08	14/23	2	606.00	1212
FBG6319AS	200	0.10	14/23	2	849.00	1698
FBG6319BS	250	0.12	14/23	2	854.00	1708
FBG6319CS	500	0.20	19/26	2	994.00	1988
FBG6319ES	1000	0.30	24/29	1	1455.00	1455

* Flask in Test Tube form


NEW

Flasks, Volumetric, Class 'A', with Stopper, Made of 3.3 Borosilicate Glass, - USP, in compliance with ASTM (Trapezoidal Shape)



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG65194S	5	0.02	10/19	2	949.50	1899
FBG65195S	10	0.02	10/19	2	809.00	1618


Flasks, Volumetric, Class 'A', Amber, with Penny Head Stopper, Made of 3.3 Borosilicate Glass, - USP, in compliance with ASTM



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63094S	5	0.02	10/19	2	811.00	1622
FBG63095S	10	0.02	10/19	2	704.00	1408
FBG63092S	20	0.03	10/19	2	849.00	1698
FBG63096S	25	0.03	10/19	2	781.00	1562
FBG63097S	50	0.05	10/19	2	752.00	1504
FBG63098S	100	0.08	14/23	2	733.00	1466
FBG6309AS	200	0.10	14/23	2	1033.00	2066
FBG6309BS	250	0.12	14/23	2	1038.00	2076
FBG6309CS	500	0.20	19/26	2	1121.00	2242
FBG6309ES	1000	0.30	24/29	1	1635.00	1635

NEW

Flasks, Volumetric, Class 'A', Amber, with Penny Head Stopper, Made of 3.3 Borosilicate Glass, - USP, in compliance with ASTM (Trapezoidal Shape)



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG64094S	5	0.02	10/19	2	1120.50	2241
FBG64095S	10	0.02	10/19	2	961.00	1922
FBG64092S	20	0.03	10/19	2	960.50	1921

Flask, Volumetric, Wide Mouth Class 'A', with Penny Head Stopper, in compliance with ISO




Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62196W	25	0.04	14/23	2	684.00	1368
FBG62197W	50	0.06	14/23	2	665.00	1330
FBG62198W	100**	0.10	19/26	2	699.00	1398
FBG6219AW	200	0.15	19/26	2	902.00	1804

** Not covered under standards



Glassware


Flask, Volumetric, Wide Mouth
Class 'A', Amber, with Penny
Head Stopper, in compliance
with ISO



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62705W	25	0.04	14/23	2	781.00	1562
FBG62707W	50	0.06	14/23	2	752.00	1504
FBG62708W	100**	0.10	19/26	2	733.00	1466
FBG6270AW	200	0.15	19/26	2	1033.00	2066

** Not covered under standards

Flasks, Volumetric, Class 'A',
clear, with Penny Head Stopper,
in compliance with ISO



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62190S	1*	0.015	10/19	2	789.50	1579
FBG62191S	2*	0.015	10/19	2	805.00	1610
FBG62193S	3	0.015	10/19	2	789.50	1579
FBG62194S	5	0.025	10/19	2	723.00	1446
FBG62195S	10	0.025	10/19	2	626.00	1252
FBG62192S	20	0.04	10/19	2	626.00	1252
FBG62196S	25	0.04	10/19	2	577.00	1154
FBG6219H	40**	0.06	10/19	2	634.00	1268
FBG62197S	50	0.06	12/21	2	553.00	1106
FBG62199	75**	0.06	14/23	2	756.50	1513
FBG62198S	100	0.10	14/23	2	538.00	1076
FBG6219D	150	0.10	14/23	2	1171.00	2342
FBG6219AS	200	0.15	14/23	2	757.00	1514
FBG6219BS	250	0.15	14/23	2	761.00	1522
FBG6219CS	500	0.25	19/26	2	994.00	1988
FBG6219G	900**	0.10	19/26	1	1586.00	1586
FBG6219ES	1000	0.40	24/29	1	1455.00	1455
FBG6219FS	2000	0.60	29/32	1	2363.00	2363

* Flask in Test Tube form


** Not covered under standards

Flasks, Volumetric, Class 'A',
clear, with Penny Head Stopper, in
compliance with ISO (Trapezoidal Shape)




Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG64194S	5	0.025	10/19	2	854.00	1708
FBG64195S	10	0.025	10/19	2	726.00	1452
FBG64192S	20	0.04	10/19	2	726.00	1452

Flasks, Volumetric, Class 'B',
clear, with Penny Head Stopper,
in compliance with ISO




Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62204S	5	0.05	10/19	2	223.00	446
FBG62205S	10	0.05	10/19	2	213.00	426
FBG62202S	20	0.08	10/19	2	213.00	426
FBG62206S	25	0.08	10/19	2	209.00	418
FBG62207S	50	0.12	12/21	2	228.00	456
FBG62208S	100	0.20	14/23	2	257.00	514
FBG6220AS	200	0.30	14/23	2	296.00	592
FBG6220BS	250	0.30	14/23	2	305.00	610
FBG6220CS	500	0.50	19/26	2	407.00	814
FBG6220ES	1000	0.80	24/29	1	611.00	611
FBG6220FS	2000	1.20	29/32	1	1155.00	1155

Flasks, Volumetric, Class 'B',
clear, with Penny Head Stopper,
in compliance with ISO (Trapezoidal Shape)



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG63204S	5	0.05	10/19	2	269.50	539
FBG63205S	10	0.05	10/19	2	244.00	488
FBG63202S	20	0.08	10/19	2	243.50	487

Flask, Volumetric, Class 'A',
Amber, with Penny Head
Stopper, in compliance with ISO



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/ Case	Price in ₹. per Pc	Price in ₹. per case
FBG62704S	5	0.025	10/19	2	864.00	1728
FBG62705S	10	0.025	10/19	2	752.00	1504
FBG62702S	20	0.04	10/19	2	752.00	1504
FBG62706S	25	0.04	10/19	2	694.00	1388
FBG6219I	40**	0.06	10/19	2	729.00	1458
FBG62707S	50	0.06	12/21	2	660.00	1320
FBG6219J	75**	0.06	12/21	2	1186.50	2373
FBG62708S	100	0.10	14/23	2	650.00	1300
FBG6219K	150	0.10	14/23	2	1529.00	3058
FBG6270AS	200	0.15	14/23	2	907.00	1814
FBG6270BS	250	0.15	14/23	2	888.00	1776
FBG6270CS	500	0.25	19/26	2	1194.00	2388
FBG6270ES	1000	0.40	24/29	1	1742.00	1742
FBG6270FS	2000	0.60	29/32	1	2833.00	2833

** Not covered under standards


NEW

Flask, Volumetric, Class 'A', Amber, with Penny Head Stopper, in compliance with ISO (Trapezoidal Shape)



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63704S	5	0.025	10/19	2	980.00	1960
FBG63705S	10	0.025	10/19	2	865.50	1731
FBG63702S	20	0.04	10/19	2	858.50	1717


Flask, Volumetric, Class 'B', Amber, with Penny Head Stopper, in compliance with ISO



Part No.	Capacity (ml)	Tolerance (±mm)	Stopper size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62714S	5	0.05	10/19	2	339.00	678
FBG62715S	10	0.05	10/19	2	320.00	640
FBG62712S	20	0.05	10/19	2	320.00	640
FBG62716S	25	0.08	10/19	2	305.00	610
FBG62717S	50	0.12	12/21	2	311.00	622
FBG62718S	100	0.20	14/23	2	344.00	688
FBG6271AS	200	0.30	14/23	2	398.00	796
FBG6271BS	250	0.30	14/23	2	407.00	814
FBG6271CS	500	0.50	19/26	2	582.00	1164
FBG6271ES	1000	0.80	24/29	1	922.00	922
FBG6271FS	2000	0.80	29/32	1	1567.00	1567

NEW

Flask, Volumetric, Class 'B', Amber, with Penny Head Stopper, in compliance with ISO (Trapezoidal Shape)



Part No.	Capacity (ml)	Tolerance (± ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63714S	5	0.05	10/19	2	407.50	815
FBG63715S	10	0.05	10/19	2	376.00	752
FBG63712S	20	0.05	10/19	2	375.50	751

Funnels, Plain, 60°C Angle, Long Stem



Part No.	Diameter (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62221	25	10	83.60	836
FBG62222	50	10	100.40	1004
FBG62223	65	10	122.80	1228
FBG62224	75	10	122.00	1220
FBG62225	100	10	166.90	1669
FBG62226	125	2	335.50	671

Funnels, Separating, Straightbore, Glass Key Stopcock, Pear Shape



Part No.	Capacity (ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62237	60	14/19	1	764.00	764
FBG62238	100	19/26	1	758.00	758
FBG6223B	250	29/32	1	970.00	970
FBG6223C	500	29/32	1	1304.00	1304
FBG6223E	1000	29/32	1	1635.00	1635

Funnels, Separating, Rotaflow Stopcock, Pear Shape, PTFE Key



Part No.	Capacity (ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62247	60	14/19	1	764.00	764
FBG62248	100	19/26	1	758.00	758
FBG6224B	250	29/32	1	970.00	970
FBG6224C	500	29/32	1	1304.00	1304
FBG6224E	1000	29/32	1	1635.00	1635

Funnels, Separating, Straight Bore, PTFE Stopcock, Pear Shape



Part No.	Capacity (ml)	Stopper Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62257	60	14/19	1	792.00	792
FBG62258	100	19/26	1	882.00	882
FBG6225B	250	29/32	1	1215.00	1215
FBG6225C	500	29/32	1	1308.00	1308
FBG6225E	1000	29/32	1	1568.00	1568

Gas Generator, Kips. With I/C Joints




Part No.	Capacity (ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6227C	500	1	7567.00	7567
FBG6227E	1000	1	9752.00	9752




Glassware

Museum Jars



Part No.	Diameter X Height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6252A	100 x 200	2	2436.00	4872
FBG6252B	100 x 250	2	2690.50	5381
FBG6252C	80 x 300	2	2834.50	5669
FBG6252D	100 x 300	2	3268.50	6537
FBG6252E	100 x 375	2	3481.50	6963
FBG6252F	80 x 380	2	3923.50	7847
FBG6252G	150 x 300	2	4713.50	9427
FBG6252H	150 x 380	2	5385.00	10770
FBG6252J	220 x 300	2	11026.00	22052

Pipette, Milk Testing, Gerber, in compliance with ISO




Part No.	Capacity (ml)	Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62368	10.75	140	10	140.30	1403

Petridishes



Part No.	Diameter x height (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6329A	90 x 15	72	133.17	9588
FBG6329B	100 x 15	72	153.69	11066
FBG6329C	100 x 20	72	174.42	12558

Pipettes, Graduated, Serological, Class 'A', Colour Coded, For Complete Delivery, in compliance with ISO




Part No.	Capacity (ml)	Sub division (ml)	Tolerance (± ml)	Colour	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63301	1	0.01	0.006	Yellow	10	348.00	3480
FBG63302	2	0.02	0.010	Black	10	325.00	3250
FBG63303	5	0.05	0.030	Red	10	155.00	1550
FBG63304	10	0.10	0.050	Orange	10	455.80	4558
FBG63306	25	0.20	0.100	White	10	534.70	5347

Pipettes For Milk Testing, in compliance with ISO




Part No.	Capacity (ml)	toler-ance (± ml)	Marked at (ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62361	1.1	0.025	0.5, 1.1	10	140.20	1402
FBG62362	2.2	0.040	1.0, 2.0, 2.1 & 2.2	10	151.50	1515

Pipettes, Graduated, Serological, Class 'B', Colour Coded, For Complete Delivery, in compliance with ISO



Part No.	Capacity (ml)	Sub division (ml)	Tolerance (± ml)	Colour	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63311	1	0.01	0.01	Yellow	20	111.70	2234
FBG63312	2	0.02	0.02	Black	20	128.70	2574
FBG63313	5	0.05	0.05	Red	20	140.10	2802
FBG63314	10	0.10	0.10	Orange	20	165.00	3300
FBG63316	25	0.20	0.20	White	20	224.40	4488

Pipettes, Graduated Serological, Class 'B', White colour Graduation, in compliance with ISO




Part No.	Capacity (ml)	Sub division (ml)	Tolerance (\pm ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62311	1	0.01	0.01	10	117.40	1174
FBG62312	2	0.02	0.02	10	123.20	1232
FBG62313	5	0.05	0.05	10	134.50	1345
FBG62314	10	0.10	0.10	10	145.80	1458

Pipettes Pasteur




Part No.	Size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62363	144	100	14.54	1454
FBG62364	240	100	16.81	1681

Pipette, Measuring, Mohr Type, Amber Graduation, Class 'A', in compliance with ISO




Part No.	Capacity (ml)	Sub division (ml)	Tolerance (\pm ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63341	1	0.01	0.006	10	348.00	3480
FBG63342	2	0.02	0.010	10	359.30	3593
FBG63343	5	0.05	0.030	10	427.50	4275
FBG63344	10	0.10	0.050	10	450.10	4501
FBG63346	25	0.20	0.100	10	529.00	5290

Pipettes, Transfer, Volumetric, with Colour Code, Class 'A', in compliance with ISO




Part No.	Capacity (ml)	Tolerance (\pm ml)	Colour Code	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63321	1	0.007	Blue	10	505.00	5050
FBG63322	2	0.010	Orange	10	519.00	5190
FBG63323	5	0.015	White	10	563.00	5630
FBG63324	10	0.020	Red	10	359.50	3595
FBG63325	20	0.030	Yellow	10	461.30	4613
FBG63326	25	0.030	Blue	10	961.00	9610
FBG63327	50	0.050	Red	10	561.50	5615
FBG63328	100	0.080	Yellow	5	788.20	3941

Pipette, Measuring, Mohr Type, Class 'B', in compliance with ISO



Part No.	Capacity (ml)	Sub division (ml)	Tolerance (\pm ml)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6235A	0.1	0.01	0.01	10	123.20	1232
FBG6235B	0.2	0.01	0.01	10	123.20	1232
FBG62351	1	0.01	0.01	10	117.40	1174
FBG62352	2	0.02	0.02	10	123.20	1232
FBG62353	5	0.05	0.05	10	134.50	1345
FBG62354	10	0.10	0.10	10	145.80	1458
FBG62356	25	0.20	0.20	10	224.40	2244

Pipettes, Transfer, Volumetric, with Colour Code, Class 'A', in compliance with ASTM (New)



Part No.	Capacity (ml)	Tolerance (\pm ml)	Colour Code	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG64321	1	0.006	Blue	10	260.00	2600
FBG64322	2	0.006	Orange	10	325.00	3250
FBG64329	3	0.010	Yellow	10	339.00	3390
FBG64330	4	0.010	Orange	10	297.40	2974
FBG64323	5	0.010	White	10	344.00	3440
FBG64331	6	0.010	Blue	10	340.50	3405
FBG64332	7	0.010	Orange	10	340.30	3403
FBG64333	8	0.020	Red	10	340.30	3403
FBG64334	9	0.020	White	10	340.30	3403
FBG64324	10	0.020	Red	10	351.60	3516
FBG64324A	15	0.030	Orange	10	394.70	3947
FBG64325	20	0.030	Yellow	10	458.90	4589
FBG64326	25	0.030	Blue	10	491.30	4913
FBG64326A	30	0.030	Yellow	10	490.80	4908
FBG64326B	40	0.050	White	10	512.20	5122
FBG64327	50	0.050	Red	10	553.30	5533
FBG64328	100	0.080	Yellow	10	844.40	8444



Glassware

Pipettes, Transfer, Volumetric, with Colour Code, Class 'B', in compliance with ISO



Part No.	Capacity (ml)	Tolerance (± ml)	Colour Code	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62331	1	0.015	Blue	10	117.40	1174
FBG62332	2	0.020	Orange	10	117.40	1174
FBG62333	5	0.030	White	10	123.20	1232
FBG62334	10	0.040	Red	10	145.80	1458
FBG62335	20	0.060	Yellow	10	173.80	1738
FBG62336	25	0.060	Blue	10	179.50	1795
FBG62337	50	0.100	Red	10	201.40	2014
FBG62338	100	0.160	Yellow	5	324.20	1621

Tubes, Centrifuge, Graduated, with I/C Stopper



Part No.	Capacity (ml)	App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6245A	10	15 x 110	25	304.56	7614
FBG62458	15	17 x 120	25	321.56	8039
FBG62459	50	29 x 135	10	388.30	3883

Tubes, Centrifuge, Conical Bottom, Plain



Part No.	Capacity (ml)	Dimension App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62451	15	17 x 120	100	46.93	4693

Tubes, Flat Bottom, with Screw Cap



Part No.	Capacity (ml)	Dimension App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62463	5	18 x 45	100	38.00	3800
FBG62464	15	25 x 57	100	49.00	4900
FBG62466	30	25 x 100	100	63.00	6300
FBG62467	60	25 x 200	100	89.14	8914

Tubes, Centrifuge, Plain, with I/C Stopper



Part No.	Capacity (ml)	App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62452	15	17 x 120	25	304.88	7622

Tubes, Centrifuge, Graduated, Conical Bottom, Amber Stain Graduation



Part No.	Capacity (ml)	Dimension App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62455	15	17 x 120	25	116.00	2900
FBG62457	50	29 x 135	25	212.60	5315

Tubes, Amber, Flat Bottom, with Screw Cap



Part No.	Capacity (ml)	Dimension App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63463	5	18 x 45	100	83.90	8390
FBG63464	15	25 x 57	100	116.00	11600
FBG63466	30	25 x 100	100	145.00	14500
FBG63467	60	25 x 200	100	144.97	14497

Tubes, Round Bottom, with Screw Cap



Part No.	Capacity (ml)	Dimension App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG62493	5	18 x 45	100	44.60	4460
FBG62494	15	25 x 57	100	44.41	4441
FBG62496	30	25 x 100	100	63.00	6300
FBG62497	60	25 x 200	100	77.79	7779
FBG62499	150	38 x 200	50	286.00	14300

Test Tubes, Graduated, with I/C Stopper




Part No.	Capacity (ml)	O.D. x Length (mm)	Stopper size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6251C	5	12 x 100	10/15	25	139.92	3498
FBG6251D	10	15 x 125	12/14	25	185.08	4627
FBG6251H	25	22 x 150	14/15	25	150.00	3750
FBG6251J	30	22 x 175	14/15	25	207.04	5176
FBG6251K	50	25 x 200	19/20	10	228.00	2280
FBG6251L	100	32 x 200	24/25	10	364.10	3641

Tubes, Amber, Round Bottom, with Screw Cap



Part No.	Capacity (ml)	Dimension App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG63493	5	18 x 45	100	78.23	7823
FBG63494	15	25 x 57	100	106.11	10611
FBG63496	30	25 x 100	100	128.56	12856
FBG63497	60	25 x 200	100	156.32	15632
FBG63498	150	38 x 200	50	659.34	32967

Test Tube, Borosilicate, w/o Rim




Part No.	App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6347A	10 x 75	100	8.48	848
FBG6347B	12 x 75	100	8.56	856
FBG6347C	12 x 100	100	10.92	1092
FBG6347D	15 x 125	100	13.33	1333
FBG6347E	15 x 150	100	15.77	1577
FBG6347F	16 x 100	100	14.55	1455
FBG6347G	18 x 150	100	18.17	1817
FBG6347H	25 x 100	100	23.97	2397
FBG6347J	25 x 150	100	26.68	2668
FBG6347K	25 x 200	100	38.35	3835
FBG6347L	32 x 200	100	60.26	6026
FBG6347M	38 x 200	100	78.53	7853

Tubes, Test, with I/C Stopper



Part No.	Capacity (ml)	O.D x Length (mm)	Stopper size (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6250C	5	12 x 100	10/15	25	106.24	2656
FBG6250D	10	15 x 125	12/14	25	123.04	3076
FBG6250H	25	22 x 150	14/15	25	156.44	3911
FBG6250K	50	25 x 200	19/20	10	195.20	1952
FBG6250L	100	32 x 200	24/25	10	313.50	3135

Test Tube, Borosilicate, with Rim



Part No.	App. O.D x Length (mm)	Qty/Case	Price in ₹. per Pc	Price in ₹. per case
FBG6348A	10 x 75	100	8.48	848
FBG6348B	12 x 75	100	8.56	856
FBG6348C	12 x 100	100	10.92	1092
FBG6348D	15 x 125	100	13.33	1333
FBG6348E	15 x 150	100	15.77	1577
FBG6348G	18 x 150	100	18.17	1817
FBG6348H	25 x 100	100	23.10	2310
FBG6348J	25 x 150	100	25.43	2543
FBG6348K	25 x 200	100	38.35	3835
FBG6348L	32 x 200	100	60.26	6026
FBG6348M	38 x 200	100	78.53	7853



Technical Data

Low Expansion Borosilicate Glass

Form the 16th Century to, today, chemical research teams have used glass containers for a very basic reason the glass containers is transparent, almost invisible. And so the contents and the reaction are clearly visible, But because chemists must heat, cool and mix chemical substances, ordinary glass is not always adequate for laboratory works.

Laboratory works requires apparatus made in a glass - which can readily be moulded into any desired shape or form, which offers maximum inertness when in contact with the widest range of chemical substances, which can withstand thermal shock with fracture and high temperature work without deforming, and which will be resilient enough to survive the everyday knocks to which it will be subjected in normal laboratory handling, washing and sterilizing processes.

Chemical Composition

Fisher^{brand} Glassware is a low alkali borosilicate composition. Its typical chemical composition is given under. It is virtually free of magnesia-lime-Zinc group and contains only traces of heavy metals.

Thermal Properties

As the Coefficient of thermal expansion of Borosilicate glass is low, the thermal stresses under a given temperature gradient are consequently low and the glass can withstand higher temperature gradients and also sudden temperature changes/thermal shocks. Minute scratching of glass surface can however reduce its thermal resistance. In general the 'Strain Point' should be regarded as the maximum safe operating temperature of Fisher^{brand} glassware. When heated above 500°C the glass may acquire permanent stresses on cooling. All Fisher^{brand} labware is annealed in modern ovens under strictly controlled conditions to ensure minimal residual stress in the products.

Chemical Durability

Fisher^{brand} Glassware in highly resistance to water, neutral and acid solutions, concentrated acids and their mixtures as well as to chloride, bromine, iodine, and organic matters. Even during extended period of reaction and at temperatures above 100°C, its chemical resistance exceeds that of most metals and other materials. It can withstand repeated dry and wet sterilisation without surface deterioration and subsequent contamination. Resistance to attack of various chemicals is shown under. Only hydrofluoric acid, very hot phosphoric acid and alkaline solutions increasingly attack the glass surface with rising concentration and temperature.

Fabrication with Borosilicate Glass

Due to low expansion of glass and easy workability, this glass can be shaped, formed, joined into complicated apparatus. It can be done even by an analyst in his own laboratory. He can keep on changing till he gets what he needs. In case where annealing in a controlled oven is difficult he can do so by flame annealing which is also great advantage.

Optical Properties

Laboratory glassware made from Borosilicate Glass shows no noticeable absorption in the visible region of the spectrum. It appear consequently clear and colourless.

When treated with proper care Fisher^{brand} laboratory apparatus will give a long and satisfactory service. The following prepared notes are to assist users in obtaining the maximum life and performance from their apparatus. Our sales department will be happy to advise on any aspect concerning the safe use of our products.



PERCENTAGE BY WEIGHT	
SiO ₂	81%
B ₂ O ₃	13%
Na ₂ O	4%

Coefficient of Linear Expansion	32.5x10⁻⁷/°C
Strains Point	515°C
Annealing Point	565°C
Softening Point	820°C
Specific Heat	0.2
Thermal Conductivity (Cal/cm ² /°C/Sec)	0.0027

Contact Chemical	Duration In Hour	Loss in WT. MG/M ²
Water distilled at 100°C	6	10
Water Vapour Steam at 121°C	1	75
Acid HCl	6	100
80% H ₂ SO ₄ at 130°C	12	140
Alkali - 1N soln. of Na ₂ CO ₃ boiling	6	4000
Infusion fluids Isotonic		
NaCl (0.85%) 121°C	2.5	70
Glucose (5%) 121°C	2.5	50





Fisherbrand glassware assures accuracy and precision in volumetric measurements

- **Volumetric Flasks**
Clear, Amber (ISO / ASTM)
- **Pipettes**
Graduated, Volumetric, Mohr (ISO / ASTM)
- **Measuring Cylinders**
Hexagonal, Stoppered (ISO / ASTM)
- **Burettes**
Glass Key, PTFE Key (ISO)

Other Key Products in Glassware Range

- **Bottles**
Clear, Amber
- **Condensers**
Allihn, Liebig, Graham
- **Desiccators**
Vacuum, Porcelain Plate
- **Tubes**
Clear, Amber
Rim, Rimless
- **Beakers**
Low Form
Tall Form
- **Funnels**
- **Weighing Scoop**
- **Flasks**
Round/Flat Bottom
2 or 3 necked



Easy tests for all purposes



Test papers and test strips



Colorimetric and titrimetric test kits



Filter oaoers and extraction thimbles

MACHEREY-NAGEL
www.mn-net.com



Cat. No.	SAP Code	Product Name	Pack size	Dimensions	Price in ₹
Filter Papers					
N1120	10541943	Filter Aids-2101	500 GM	500gm	10476
88011	10221921	615A equivalent to grade 1	100 Circle	125mm	969
88151	Q88151	640W equivalent to grade 41	100 Circle	125mm	3377
88061	Q88061	617 equivalent to grade 4	100 Circle	125mm	1265
88071	Q88071	618 equivalent to grade 3	100 Circle	110mm	1398
88021	10025524	615A equivalent to grade 1	100 Sheets	46X57 cm	9195
88155	10511551	640md equivalent to grade 40	100 Circle	120mm	3074
88001	12553627	615A equivalent to grade 1	100 Circle	110mm	804
88195	10024240	640de equivalent to grade 42	100 Circle	125mm	3184
89391	10032691	1670 equivalent to Grade 54	100 Circle	110mm	3468
88081	10086491	618 equivalent to grade 3	100 Circle	125mm	2018
88031	12533657	615A equivalent to grade 1	100 Sheets	58X58 cm	12655
89414	Q89414	1672 equivalent to Grade 52	100 Circle	110mm	3986
89415	10073950	1672 equivalent to Grade 52	100 Circle	125mm	4429
88212	11512391	1640m equivalent to Grade 541	100 Circle	150mm	5573
88002	12563627	615A equivalent to grade 1	100 Circle	90mm	790
89444	Q89444	1674 equivalent to Grade 50	100 Circle	110mm	3751
88012	12553637	615A equivalent to grade 1	100 Circle	150mm	1217
89450	Q89450	1674 equivalent to Grade 50	100 Circle	320mm	21897
88216	11518040	1640d equivalent to Grade 540	100 Circle	150mm	5543
88193	10544031	640de equivalent to grade 42	100 Circle	90mm	2123
89445	Q89445	1674 equivalent to Grade 50	100 Circle	125mm	3995
88211	11502391	1640m equivalent to Grade 541	100 Circle	125mm	3500
88152	Q88152	640w equivalent to grade 41	100 Circle	150mm	4338
88153	10712331	640md equivalent to grade 40	100 Circle	90mm	1941
88154	10046471	640md equivalent to grade 40	100 Circle	110mm	2524
88156	10442491	640md equivalent to grade 40	100 Circle	150mm	4085
88121	10071712	619de equivalent to grade 5	100 Circle	125mm	1534
88196	Q88196	640de equivalent to grade 42	100 Circle	150mm	4269
88140	10261531	640w equivalent to grade 41	100 Circle	90mm	2239
86031	10013151	Chromatography paper MN261 equivalent to Grade 20 CHR	25 Sheets	58X60 cm	15274
88194	10460381	640de equivalent to grade 42	100 Circle	110mm	2711
88051	12593667	617 equivalent to grade 4	100 Circle	110mm	1243
88695	Q88695	619eh equivalent to grade 6	100 Circle	125mm	1767



Filter Papers


Cat. No.	SAP Code	Product Name	Pack size	Dimensions	Price in ₹
89455	89455	1674 equivalent to Grade 50	100 Circle	58X58 cm	31248
89395	89395	1670 equivalent to Grade 54	100 Circle	58X58 cm	61705
88215	11567920	1640 d equivalent to Grade 540	100 Circle	125mm	5201
88235	Q88235	1640de equivalent to Grade 542	100 Circle	125mm	5135
88236	88236	1640de equivalent to Grade 542	100 Circle	150mm	5640
729227	729227	Chromafil Xtra CA-45/25	100/ Pack	25 mm	41896
M203009	M203009	640m equivalent to grade 43	100 Circle	90mm	2651
M203011	M203011	640m equivalent to grade 43	100 Circle	110mm	3384
M203012	M203012	640m equivalent to grade 43	100 Circle	125mm	3615
M203015	M203015	640m equivalent to grade 43	100 Circle	150mm	5399
M206009	M206009	640dd equivalent to grade 44	100 Circle	90mm	2535
M206011	M206011	640dd equivalent to grade 44	100 Circle	110mm	3836
M206012	M206012	640dd equivalent to grade 44	100 Circle	125mm	4498
M206015	M206015	640dd equivalent to grade 44	100 Circle	150mm	5631
88141	Q88141	640w equivalent to grade 41	100 Circle	110mm	2893
88065	11342711	616md equivalent to grade 2	100 Circle	125mm	1175
203012	M203012	Filter Paper Circles MN 640 m 12,5 cm diameter OEM presentation pack of 100	100 Circle	125mm	3615

Tips to Select the Right Filter Paper

To ensure effective filtration the choice depends on:

- The filtration method to be used.
- The filter paper grade.
- The properties of the medium to be filtered.
- The size of the particulate solids to be removed
- The required degree of clarification

Cat. No.	SAP Code	Product Name	Range	Pack size	Price in ₹
ECO					
M1008	M1008	Ammonium	(0.2-3.0 Mg/L)	50 Tests	9450
M1208	M1208	Ammonium -Refill	(0.2-3.0 Mg/L)	50 Tests	4725
M1018	10085701	Chloride	(1.0-60 Mg/L)	90 Tests	9450
M1218	M1218	Chloride-Refill	(1.0-60 Mg/L)	90 Tests	6825
M1015	11419578	Chlorine	(0.1-2 Mg/L)	150 Tests	6668
M1215	M1215	Chlorine-Refill	(0.1-2 Mg/L)	150 Tests	3990
M1020	10157220	Chromium (Vi)	(0.02-0.5 Mg/L)	140 Tests	8190
M1220	10366212	Chromium (Vi)-Refill	(0.02-0.5 Mg/L)	140 Tests	6229
M1037	11419588	Copper	(0.1-1.5 Mg/L)	100 Tests	7875
M1237	M1237	Copper-Refill	(0.1-1.5 Mg/L)	100 Tests	4935
M1022	11459578	Cyanide	(0.01-0.20 Mg/L)	100 Tests	12238
M1023	M1023	Cyanuric Acid	(10-100 Mg/L)	100 Tests	12863
M1223	M1223	Cyanuric Acid-Refill	(10-100 Mg/L)	100 Tests	6123
M1026	11489578	Iron	(0.04-1 Mg/L)	100 Tests	8400
M1226	10769224	Iron-Refill	(0.04-1 Mg/L)	100 Tests	5250
M1038	10708654	Manganese	(0.1-1.5 Mg/L)	70 Tests	10658
M1238	M1238	Manganese-Refill	(0.1-1.5 Mg/L)	70 Tests	5959
M1040	10457292	Nickel	(0.1-1.5 Mg/L)	150 Tests	13064
M1240	13448205	Nickel-Refill	(0.1-1.5 Mg/L)	150 Tests	6930
M1041	10008910	Nitrate	(4-120 Mg/L)	120 Tests	6615
M1241	M1241	Nitrate-Refill	(4-120 Mg/L)	120 Tests	3528
M1044	11429588	Nitrite	(0.02-0.5 Mg/L)	150 Tests	7552
M1244	M1244	Nitrite-Refill	(0.02-0.5 Mg/L)	150 Tests	3833
M1088	M1088	Oxygen	(1-10 Mg/L)	50 Tests	8715
M1288	M1288	Oxygen-Refill	(1-10 Mg/L)	50 Tests	6090
M1084	11499588	Phosphate	(0.2-5 Mg/L)	100 Tests	8190
M1284	M1284	Phosphate-Refill	(0.2-5 Mg/L)	100 Tests	4673
M1033	10171852	Silica / Silicon	(0.2-3.0 Mg/L)	100 Tests	8902
M1233	AAB-M1233	Silica- / Silicon-Refill	(0.2-3.0 Mg/L)	100 Tests	5460
HE					
M0006	10015220	Ammonium	(0.02-0.5 Mg/L)	110 Tests	15750
M0106	M0106	Ammonium -Refill	(0.02-0.5 Mg/L)	110 Tests	8925
M0015	M0015	Chlorine(Free + Total)	(0.02-0.60 Mg/L)	160 Tests	14280
M0115	M0115	Chlorine(Free + Total)-Refill	(0.02-0.60 Mg/L)	160 Tests	9426
M0050	10470761	Copper	(0.04-0.5 Mg/L)	150 Tests	15750
M0150	M0150	Copper-Refill	(0.04-0.5 Mg/L)	150 Tests	9426
M0028	10075240	Cyanide	(0.002-0.04 Mg/L)	55 Tests	14280
M0128	M0128	Cyanide-Refill	(0.002-0.04 Mg/L)	55 Tests	8820
M0040	10024260	Iron	(0.01-0.2 Mg/L)	300 Tests	14897
M0140	M0140	Iron-Refill	(0.01-0.2 Mg/L)	300 Tests	8925
M0055	M0055	Manganese	(0.03-0.5 Mg/L)	100 Tests	14018
M0155	M0155	Manganese-Refill	(0.03-0.5 Mg/L)	100 Tests	8820
M0063	M0063	Nitrite	(0.005-0.1 Mg/L)	150 Tests	16065
M0163	M0163	Nitrite-Refill	(0.005-0.1 Mg/L)	150 Tests	9345
M0082	M0082	Phosphate	(0.05-1 Mg/L)	300 Tests	14280
M0182	M0182	Phosphate-Refill	(0.05-1 Mg/L)	300 Tests	8925
M0080	M0080	Phosphate (Dev)	(0.01-0.25 Mg/L)	100 Tests	14018
M0180	M0180	Phosphate (Dev)-Refill	(0.01-0.25 Mg/L)	100 Tests	8820
M0087	10218180	Silica / Silicon	(0.01-0.3 Mg/L)	120 Tests	14280
M0187	10641363	Silica / Silicon-Refill	(0.01-0.3 Mg/L)	120 Tests	9030




Cat. No.	SAP Code	Product Name	Range	Pack size	Price in ₹
PEHANON INDICATOR PAPERS					
M0401	10154623	Ph 1.0-12.0	(Gradations 1 -2 -3 -5 -6 -7 -8 -9 -10 -12)	200 Tests	3208
M0411	M0411	Ph 0.0-1.8	(Gradations 0-0.3-0.6 -0.8 -1.0 -1.2 -1.4 -1.6 -1.8)	200 Tests	4331
M0412	10460431	Ph 1.0-2.8	(Gradations 1.0 -1.3 -1.6 -1.8 -2.0 -2.2 -2.5 -2.8)	200 Tests	4331
M0413	10330961	Ph 1.8-3.8	(Gradations 1.8 -2.1 -2.4 -2.7 -3.0 -3.32 -3.5 -3.8)	200 Tests	2310
M0414	M0414	Ph 2.8-4.6	(Gradations 2.8 -3.1 -3.4 -3.6 -3.8 -4.0 -4.3 -4.6)	200 Tests	2370
M0415	10100061	Ph 3.8-5.5	(Gradations 3.8 -4.0 -4.2 -4.4 -4.6 -4.9 -5.2 -5.5)	200 Tests	2370
M0424	10547502	Ph 4.0-9.0	(Gradations 4.0 -4.5 -5.0 -5.5 -6.0 -6.5 -7.0 -7.5 -8.0 -8.5 -9.0)	200 Tests	3150
M0416	10196160	Ph 5.2-6.8	(Gradations 5.2 -5.5 -5.7 -5.9 -6.1 -6.3 -6.5 -6.8)	200 Tests	2480
M0417	10708911	Ph 6.0-8.1	(Gradations 6.0 -6.3 -6.6 -6.9 -7.2 -7.5 -7.8 -8.1)	200 Tests	2415
M0419	10689572	Ph 7.2-8.8	(Gradations 7.2 -7.4 -7.6 -7.8 -8.0 -8.2 -8.5 -8.8)	200 Tests	2480
M0420	10736991	Ph 8.0-9.7	(Gradations 8.0 -8.2 -8.4 -8.6 -8.8 -9.1 -9.4 -9.7)	200 Tests	2480
M0421	10687901	Ph 9.5-12.0	(Gradations 9.5 -10.0 -10.5 -11.0 -11.5 -12.0)	200 Tests	2480
M0422	M0422	Ph 10.5-13.0	(Gradations 10.5 -11.0 -11.5 -12.0 -12.5 -13.0)	200 Tests	2480
M0423	10771601	Ph 12.0-14.0	(Gradations 12.0 -12.5 -13.0 -13.5 -14.0)	200 Tests	2480
M2110	11418019	Ph 0 - 14	(0 -1 -2 -3 -4 -5 -6 -7 -8 -9 -10 -11 -12 -13 -14)	100 Tests	3570
M2115	11468019	Ph 0 - 6	(0-0.5-1.0-1.5 -2.0-2.5 -3.0-3.5 -4.0 -4.5 -5.0-5.5 -6.0)	100 Tests	2429
M2118	10281701	Ph 2 - 9	(2.0-2.5-3.0-3.5-4.0-4.5-5.0-5.5-6.0-6.5-7.0-7.5-8.0-8.5-9.0)	100 Tests	3465
M2120	11488019	Ph 4.5 - 10	(4.5 -5.0 -5.5 -6.0 -6.5 -7.0 -7.5 -8.0 -8.5 -9.0 -9.5 -10.0)	100 Tests	2429
M2125	10271751	Ph 7 - 14	(7.0-7.5-8.0-8.5-9.0-9.5-10.0-10.5-11.0-11.5-12.0-12.5-13.0-13.5-14.0)	100 Tests	3438
M2180	M2180	Ph 0.3 - 2.3	(0.3 -0.7 -1.0 -1.3 -1.6 -1.9 -2.3)	100 Tests	2520
M2190	M2190	Ph 1.7 - 3.8	(1.7 -2.0 -2.3 -2.6 -2.9 -3.2 -3.5 -3.8)	100 Tests	2429
M2130	10353641	Ph 3.6 - 6.1	(3.6 -4.1 -4.4 -4.7 -5.0 -5.3 -5.6-6.1)	100 Tests	2429
M2140	M2140	Ph 5.1 - 7.2	(5.1 -5.4 -5.7 -6.0 -6.6 -6.9 -7.2)	100 Tests	2429
M2150	11498029	Ph 6.0 - 7.7	(6.0 -6.4 -6.7 -7.0 -7.3 -7.7)	100 Tests	2429
M2160	11468039	Ph 7.5 - 9.5	(7.5 -7.9 -8.2 -8.4 -8.6 -8.8 -9.1 -9.5)	100 Tests	2429
M2170	10727181	Ph 7.9 - 9.8	(7.9 -8.3 -8.6 -8.9 -9.1 -9.4 -9.8)	100 Tests	2468
Quantofix (Water testing Strips)					
M1307	10186080	Aluminium*	0 - 5 - 20 - 50 - 200 - 500	100 Tests	6947
M1315	12753884	Ammonium*	0 - 10 - 25 - 50 - 100 - 200 - 400	100 Tests	6828
M1332	11316594	Arsenic*	0 - 0.05 - 0.1 - 0.5 - 1.0 - 1.7 - 3.0	100 Tests	10768
M1334	10581363	Arsenic (Low Range)*	0-0.01 - 0.025 - 0.05 - 0.1 - 0.5100	100 Tests	12131
M1314	M1314	Ascorbic Acid	0-50-100-200-300-500-1000-2000	100 Tests	4435
M1324	10685811	Calcium*	0 - 10 - 25 - 50 - 100	100 Tests	8623
M1323	11413699	Carbonate Hardness	0 - 3 - 6 - 10 - 15 - 20 0D**	100 Tests	4049
M1321	10167430	Chloride	0 - 500 - 1000 - 1500 - 2000 - >3000	100 Tests	4894
M1317	10258390	Chlorine*	0 - 1 - 3 - 10 - 30 - 100	100 Tests	8080
M1301	10100441	Chromate*	0 - 3 - 10 - 30 - 100	100 Tests	5175
M1303	AAB-M1303	Cobalt	0 - 10 - 25 - 50 - 100 - 250 - 500 - 1000	100 Tests	4752
M1304	10536732	Copper	0 - 10 - 30 - 100 - 300	100 Tests	4536
M1318	10711811	Cyanide*	0 - 1 - 3 - 10 - 30	100 Tests	8025
M1328	11453699	Formaldehyde*	0 - 10 - 20 - 40 - 60 - 100 - 200	100 Tests	7045
M1325	11423699	Molybdenum*	0 - 5 - 20 - 50 - 100 - 250	100 Tests	9024
M1305	10100491	Nickel	0 - 10 - 25 - 50 - 100 - 250 - 500 - 1000	100 Tests	4099
M1313	10384071	Nitrate/Nitrite	0 - 10 - 25 - 50 - 100 - 250 - 500	100 Tests	4223
M1311	10474841	Nitrite	0 - 1 - 5 - 10 - 20 - 40 - 80	100 Tests	4038
M1322	10127480	Nitrite	3000 0 - 0.1 - 0.3 - 0.6 - 1.0 - 2.0 - 3.0	100 Tests	4014
M1333	10233952	Peroxide	10000 - 50 - 150 - 300 - 500 - 800 - 1000	100 Tests	4388
M1312	11493689	Peroxide	1000 - 1 - 3 - 10 - 30 - 100	100 Tests	3845
M1319	10528482	Peroxide	250 - 0.5 - 2 - 5 - 10 - 25	100 Tests	3845
M1320	10686351	Phosphate*	0 - 3 - 10 - 25 - 50 - 100	100 Tests	6120
M1316	10696351	Potassium*	0 - 200 - 400 - 700 - 1000 - 1500	100 Tests	7333

**SPECIAL INDICATOR PAPERS / UNIVERSAL INDICATOR PAPERS
VISOCOLOR KITS TITRATION**

Water Testing Kits

Cat. No.	SAP Code	Product Name	Range	Pack size	Price in ₹
M1329	11463699	Sulphate	0 - 400 - 800 - 1200 - 1600	100 Tests	4266
M1306	11463689	Sulphite	0 - 10 - 25 - 50 - 100 - 250 - 500 - 1000	100 Tests	4099
M1309	M1309	Tin	0 - 10 - 25 - 50 - 100 - 250 - 500	100 Tests	5917
M1310	11483689	Zinc*	0 - 2 - 5 - 10 - 25 - 50 - 100	100 Tests	6749
M1326	11433699	Multistick For Aquarium Total Hardness / Carbonate / Ph		100 Tests	6617
M1201	10632381	Aquadur-Water Hardness	0-5-10-15- 20 - 25 0D**	100 Tests	4102
M0734	10616371	Fluoride (Discs)*	0 - 2 - 5 - 10 - 20 - 50 - 100	30 Tests	11735
SPECIAL INDICATOR PAPERS					
M0205	10098250	Ph 0.5-5.5	(Gradations 0.5 -1.0 -1.5 -2.0 -2.5 -3.0 -3.5 -4.0 -4.5 -5.0 -5.5)	1 Reel	3032
M0225	10718881	Ph 0.5-5.5 (Refill)	(Gradations 0.5 -1.0 -1.5 -2.0 -2.5 -3.0 -3.5 -4.0 -4.5 -5.0 -5.5)	3 reels	4883
M0206	10323361	Ph 3.8-5.8	(Gradations <3.8 -3.8 -4.1 -4.3 -4.5 -4.7 -4.9 -5.2 -5.5 -5.8 ->5.8)	1 Reel	3032
M0226	M0226	Ph 3.8-5.8 (Refill)	(Gradations <3.8 -3.8 -4.1 -4.3 -4.5 -4.7 -4.9 -5.2 -5.5 -5.8 ->5.8)	3 reels	4881
M0207	10539641	Ph 4.0-7.0	(Gradations 4.0 -4.3 -4.6 -4.9 -5.2 -5.5 -5.8 -6.1 -6.4 -6.7 -7.0)	1 Reel	3032
M0227	10257840	Ph 4.0-7.0 (Refill)	(Gradations 4.0 -4.3 -4.6 -4.9 -5.2 -5.5 -5.8 -6.1 -6.4 -6.7 -7.0)	3 Reels	4881
M0208	10549641	Ph 5.4-7.0	(Gradations <5.4 -5.4 -5.7 -6.0 -6.2 -6.4 -7.0>7.0)	1 Reel	3032
M0228	M0228	Ph 5.4-7.0 (Refill)	(Gradations <5.4 -5.4 -5.7 -6.0 -6.2 -6.4 -7.0>7.0)	3 Reels	4881
M0209	10559641	Ph 5.5-9.0	(Gradations 5.5 -6.0 -6.5 -7.0 -7.5 -8.0 -8.5 -9.0)	1 Reel	3032
M0229	10241661	Ph 5.5-9.0 (Refill)	(Gradations 5.5 -6.0 -6.5 -7.0 -7.5 -8.0 -8.5 -9.0)	3 Reels	4881
M0210	10570022	Ph 6.4-8.0	(Gradations <6.4 -6.4 -6.6 -6.8 -7.0 -7.2 -7.4 -7.6 -7.8 -8.0>8.0)	1 Reel	3032
M0230	10180101	Ph 6.4-8.0 (Refill)	(Gradations <6.4 -6.4 -6.6 -6.8 -7.0 -7.2 -7.4 -7.6 -7.8 -8.0>8.0)	3 Reels	4881
M0211	10613141	Ph 7.2-9.7	(Gradations <7.2 -7.2 -7.5 -7.8 -8.1 -8.4 -8.7 -9.0 -9.3 -9.7 ->9.7)	1 Reel	3032
M0231	10404851	Ph 7.2-9.8 (Refill)	(Gradations <7.2 -7.2 -7.5 -7.8 -8.1 -8.4 -8.7 -9.0 -9.3 -9.7 ->9.7)	3 Reels	4881
M0212	10414851	Ph 8.0-10.0	(Gradations 8.0 -8.2 -8.4 -8.7 -9.0 -9.2 -9.6 -10.0)	1 Reel	3032
M0232	M0232	Ph 8.0-10.0 (Refill)	(Gradations 8.0 -8.2 -8.4 -8.7 -9.0 -9.2 -9.6 -10.0)	3 Reels	4881
M0213	10527341	Ph 9.0-13.0	(Gradations 9.0 -9.5 -10.0 -10.5 -11.0 -11.5 -12.0 -12.5 -13.0)	1 Reel	3032
M0233	10373451	Ph 9.0-13.0 (Refill)	(Gradations 9.0 -9.5 -10.0 -10.5 -11.0 -11.5 -12.0 -12.5 -13.0)	3 Reels	4881
M0214	10303651	Ph 12.0-14.0	(Gradations 12.0 -12.5 -13.0 -13.5 -14.0)	1 Reel	3032
M0234	10271851	Ph 12.0-14.0 (Refill)	(Gradations 12.0 -12.5 -13.0 -13.5 -14.0)	3 Reels	4881
UNIVERSAL INDICATOR PAPERS					
M0204	10797941	Indicator Papers Ph 1-14	Ph 1-14	1 Reel	3032
M0224	10089930	Indicator Papers Ph 1-14 (Refill)	1-14 (Refill)	3 Reels	4883
M0201	M0201	Indicator Papers Ph 1-11	Ph 1-11	1 Reel	3032
M0202	10323981	Indicator Papers Ph 1-11 (Refill)	Ph 1-11 (Refill)	3 Reels	4881
VISOCOLOR KITS TITRATION					
M5006	13215844	Titration acidity Ac 7 (Base Capacity)	(0.2-7.0 Mmol/L)	200 Tests	7508
M5206	M5206	Titration acidity Ac 7-Refill	(0.2-7.0 Mmol/L)	200 Tests	5016
M5007	10393321	Titration alkalinity Ac 7 (Total)	(0.2-7.0 Mmol/L)	200 Tests	6893
M5207	M5207	Titration alkalinity Ac 7-Refill	(0.2-7.0 Mmol/L)	200 Tests	4620
M5010	11848572	Titration calcium Ca 20	(0.1-3.6 Mmol/L)	200 Tests	8400
M5210	M5210	Titration calcium Ca 20-Refill	(0.1-3.6 Mmol/L)	200 Tests	5250
M5003	10109170	Titration carbonate Hardness C 20	(0.2-7.0 Mmol/L)	200 Tests	6893
M5203	M5203	Titration carb Hard C 20-Refill	(0.2-7.0 Mmol/L)	200 Tests	4725
M5004	10483291	Titration chloride Cl 500	(5-500 Mg/L)	300 Tests	10500
M5204	M5204	Titration chloride Cl 500-Refill	(5-500 Mg/L)	300 Tests	7875
M5009	10729994	Titration oxygen S.A 10	(0.2-10 Mg/L)	100 Tests	13127
M5209	10456861	Titration oxygen S.A 10-Refill	(0.2-10 Mg/L)	100 Tests	10011
M5008	10443671	Titration sulphite Su 100	(2-100 Mg/L)	100 Tests	11631
M5208	M5208	Titration sulphite Su 100-Refill	(2-100 Mg/L)	100 Tests	7828
M5005	12688055	Titration total Hardness H 20 F	(0.1-3.6 Mg/L)	200 Tests	7882
M5205	M5205	Titration total Hardness H 20 F-Refill	(0.1-3.6 Mg/L)	200 Tests	4601
M5002	10781283	Titration total Hardness H 2 (Residual)	(0.01-0.36 Mmol/L)	200 Tests	7882
M5202	M5202	Titration total Hardness H 2-Refill	(0.01-0.36 Mmol/L)	200 Tests	4620



Labserv® Syringe Filters / Membrane Filters

Part No.	Product Name	Pack Size	Price in ₹
LSMN2047NST	Membrane FilterNylon47 mm0.2 µmNon Sterile	100	1795
LSMN4547NST	Membrane FilterNylon47 mm0.45 µmNon Sterile	100	1795
LSMCN2047NST	Membrane FilterCellulose Nitrate47 mm0.2 µmNon Sterile	100	1901
LSMCN4547NST	Membrane FilterCellulose Nitrate47 mm0.45 µmNon Sterile	100	1901
LSMCN4547ST	Membrane FilterCellulose Nitrate47 mm0.45 µmSterile	100	2281
LSSN2013NST	Syringe FilterNylon13 mm0.2 µmNon Sterile	100	3273
LSSN4513NST	Syringe FilterNylon13 mm0.45 µmNon Sterile	100	3271
LSSN2025NST	Syringe FilterNylon25 mm0.2 µmNon Sterile	100	3377
LSSN4525NST	Syringe FilterNylon25 mm0.45 µmNon Sterile	100	3377
LSSN2013ST	Syringe FilterNylon13 mm0.2 µmSterile	100	6632
LSSN4513ST	Syringe FilterNylon13 mm0.45 µmSterile	100	6632
LSSN2025ST	Syringe FilterNylon25 mm0.2 µmSterile	50	3695
LSSN4525ST	Syringe FilterNylon25 mm0.45 µmSterile	50	3695
LSSPT2013NST	Syringe FilterPTFE13 mm0.2 µmNon Sterile	100	5752
LSSPT4513NST	Syringe FilterPTFE13 mm0.45 µmNon Sterile	100	5752
LSSPT2025NST	Syringe FilterPTFE25 mm0.2 µmNon Sterile	100	6282
LSSPT4525NST	Syringe FilterPTFE25 mm0.45 µmNon Sterile	100	6282
LSSPT2013ST	Syringe FilterPTFE13 mm0.2 µmSterile	100	7096
LSSPT4513ST	Syringe FilterPTFE13 mm0.45 µmSterile	100	7096
LSSPT2025ST	Syringe FilterPTFE25 mm0.2 µmSterile	50	4161
LSSPT4525ST	Syringe FilterPTFE25 mm0.45 µmSterile	50	4161
LSSGF1013NST	Syringe FilterGlass Fibre13 mm1.0 µmNon Sterile	100	3703
LSSGF1613NST	Syringe FilterGlass Fibre13 mm1.6 µmNon Sterile	100	3597
LSSGF1025NST	Syringe FilterGlass Fibre25 mm1.0 µmNon Sterile	100	3809
LSSGF1625NST	Syringe FilterGlass Fibre25 mm1.6 µmNon Sterile	100	3703
LSSPV2013NST	Syringe FilterPVDF13mm0.2 µmNon Sterile	100	4091
LSSPV4513NST	Syringe FilterPVDF13mm0.45 µmNon Sterile	100	4091
LSSPV2025NST	Syringe FilterPVDF25mm0.2 µmNon Sterile	100	5084
LSSPV4525NST	Syringe FilterPVDF25mm0.45 µmNon Sterile	100	5084



LabServ® Chromatography Vials

Quality performance at an affordable price

A tradition of superior quality, knowledgeable support, and innovative solutions to challenging problems make us the leader in laboratory consumables. The portfolio of Thermo Scientific™ LabServ® Autosampler Vials and Closures offers analysts and researchers:

- Patched and non-patched vials in white polypropylene boxes made from specially selected chromatographically clean materials
- Closures with different types of septa covering nearly all HPLC applications
- Convenient packs of 100 for easy laboratory use
- Combination packs for routine analytical use

Syringe & Membrane Filtration

The new LabServ® syringe-driven filter units have been developed specifically for use in sample preparation. The LabServ® syringe filters are a proven alternative for chromatography sample preparation.

LabServ® Syringe Filters

- Used for filtration of samples and solvents
- Increases the accuracy and reliability of the results
- Eliminates background peaks
- Prevents instruments' downtime
- Increases HPLC column life
- Available to suit different configurations
 - Media - Nylon, PVDF, PTFE (hydrophobic)
 - Sizes - 13mm & 25mm
 - Filtration rating - 0.2 μm , 0.45 μm



LabServ® Nylon Membrane

- Excellent compatibility with solvents
- Can be easily sterilized through steam autoclaving
- Used for filtration of microbiological media, sera & other additives
- Ideal for HPLC & Organic solvents
- Can be used for sterilizing filtration of pharmaceuticals aqueous solution & bioassay
- Used for filtration of wide range of injectable solution
- Available to suit different configurations
 - Diameter - 47mm
 - Filtration rating - 0.2 μm , 0.45 μm

LabServ® Cellulose Nitrate Membrane

- Ideal for sterility testing
- Ideal for microbial analysis
- Suitable for aqueous and protein filtration
- Used for bio-burden reduction
- Particle filtration & monitoring (high micron rating membranes)
- Available to suit different configurations
 - Diameter - 47 mm
 - Filtration rating - 0.2 μm , 0.45 μm